



**Say hello to the  
Leckey family  
of products.**

# I am delighted to introduce our new Leckey product guide and celebrate our Leckey family circle.

At Leckey we believe in **'Supporting Children'**.

To everyone in Leckey, Supporting Children means working together with our family of experts – designers, therapists, parents and the broader clinical community to achieve our goals.

Our **Mission** is; To improve the quality of life of children with special needs.

Our **Vision** is; A world where each child is included and they can access their full potential.

Our **Courageous Goal** is; To prevent hip migration in children with CP.

***In Leckey we see hip migration as the single biggest life event that our children may experience and has the greatest impact on their quality of life and their carers from that moment on.***

Everything we do in Leckey, in partnership with our family of experts, is aimed at achieving these goals.

**Our approach is simple – lead through innovation and clinical focus.**

Through early intervention, childhood right up to adulthood, Leckey's experienced team of designers and clinicians work together with our family of experts to develop exciting innovative products that make a real difference to the quality of life for the children and their families.

We aim to develop clinically and technically superior solutions, products and programs, that improve the child's Mobility, Ability and Participation.

All of our designs undergo extensive clinical and usability trials with physiotherapists, occupational therapists, users and their families.

They are supported by in-depth clinical research studies where we work with leading universities and clinicians internationally to ensure that all Leckey products deliver on the clinical needs of the healthcare professionals and the social needs of the user.

We are committed to our Mission, Vision and Courageous Goals. We want to work with you in partnership to build our family - a community of experts with shared goals and values.

**Join our family and together we will make a difference.**

**Erwin McKee**  
**General Manager Leckey, Firefly and Sunrise**  
**Medical UK Operations**



# Contents



## Seating

04

Early Activity System	06
Pal	08
Squiggles Seat	09
BeMe	10
Mygo Seat	12
Everyday Activity Seat	14
Leckey Mobility	16
Leckey Therapy Bench	18
JCM Star	20
JCM Triton	21
JCM Jupiter	22



## Standing and Walking

23

Squiggles+	24
Squiggles TT	26
Squiggles Stander	28
Mygo Stander	30
Horizon	32
Pronestander	34
Freestander	35
Totstander	36
Timoteos Ella	37
MyWay+	38
MyWay Pedal	40



## Firefly

42

GottaGo	44
GoTo Seat	46
Upsee	48
Scooot	50
Playpak	52
WeGo	54
Splashy	56
Splashy Big	58

## Other Products

59

Advance Bath Chair & Shower Trolley	60
Sleepform	61
Whitmyer Headrests	62



# Seating



# LECKEY Early Activity System

Colourful, comfortable, supportive activity system for children aged 0-36 months.



## Product Features

- Offers a gentle, playful introduction to 24-hour postural management.
- Practical all-in-one kit, providing versatile support. The positioning possibilities are endless and everything can be packaged up in the duffle bag for easy transport.
- Supplied with the popular Activity Development Programme which supports parent-led therapy at home.
- Soft-touch, anti-bacterial, machine washable and tumble dry safe covers.



## Key Positions

- Back Lying (Supine)
- Tummy Lying (Prone)
- Side Lying
- Floor Sitting (Long sitting)
- Hands and knees (Four point kneeling)



# LECKEY Pal

Colourful, comfortable classroom seating system providing mild support for children aged 1-12 years.

## Product Features

- Two-point lap belt and side pads for lateral stability. An optional ramped cushion can be added if extra help is needed to keep the pelvis and femurs in alignment.
- Where a less upright angle is required, the backrest angle and armrests adjust to suit the child's needs, and along with the large tray or table, provide appropriate support.
- Height adjustable armrests, legs and footplate.
- Available in two bright colours in wipe clean fabric [please see website for details].
- A range of base and footrest options available, including castors, rockers and stabilisers.
- Two-person and four-person matching tables available.



Available in 2 fabric colours

Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	1-3	2-5	4-8	7-12
User Height	Min 700mm [27.6"] Max 950mm [37.4"]	Min 900mm [35.4"] Max 1050mm [41.3"]	Min 1000mm [39.4"] Max 1300mm [51.2"]	Min 1200mm [47"] Max 1600mm [63"]
Max. User Weight	30kg	35kg	40kg	50kg
Seat Depth	180mm-240mm [7"-9.4"]	220mm-290mm [8.6"-11.4"]	250mm-330mm [9.8"-12.9"]	300mm-400mm [11.8"-15.7"]
Seat Width	150mm-250mm [5.9"-9.8"]	170mm-270mm [6.6"-10.6"]	230mm-330mm [9"-12.9"]	280mm-380mm [11"-14.9"]
Seat Height	160mm-240mm [6.2"-9.4"]	185mm-285mm [7.2"-11.2"]	250mm-350mm [9.8"-13.7"]	300mm-400mm [11.8"-15.7"]

# LECKEY Squiggles Seat

Suitable for use as a home chair, school chair or with a mobility base, Squiggles Seat meets all the postural and comfort requirements of a seating system for young children with special needs, while maintaining a fun, colourful design.

## Product Features

- Designed specifically for the developmental needs of younger children- the one-piece seat base and adjustable sacral support can be used to accommodate flexion or encourage extension.
- The shape of the seat and back, the laterals with their fabric tip which wrap securely around the child and the way in which the 4-point pelvic harness integrates with the hip pads, provides a snug, comfortable fit and streamlined look.
- Being low to the ground, the Squiggles Seat promotes peer to peer interaction, making it ideal for home and nursery/school environments.
- Adjustable footplate can be tilted to the preferred angle, ensuring the feet are well supported and providing a secure base for maximum upper body function. The footplate flips away to allow for easy transfers.
- Tilt-in-space feature provides a simple way to achieve a change of position without disturbing overall posture.



Available in 4 fabric colours

Technical Details	
One Size Only	
Age	1-5
Max. User Weight	22kg [48.4lbs]
Seat Depth	190mm-290mm [7.5"-11.4"]
Seat Width	160mm-260mm [6.3"-10.2"]
Chest Width	160mm-220mm [6.3"-8.75"]
Backrest Height	320mm-380mm [12.6"-15"]
Lower Leg Length	122mm-270mm [4.8"-10.6"]
Tray Height	115mm-170mm [4.5"-6.7"]
Backrest Angle	Prone 10° - Recline 25°
Seat Height (Hi-Low Chassis)	270mm-650mm [10"-25"]



# LECKEY BeMe

## Support as Individual as You Are

A revolutionary, highly configurable modular seating system which can be tailored to provide superior and individualised support, empowering a child to simply BeMe.



### Product Features

- At the heart of BeMe is the seat core, from which a range of interchangeable parts can be quickly and easily fitted. This includes three seat types (basic, contoured and split) and three backrests (basic, moderate and complex), as well as medial, lateral, proximal and distal supports, in a range of sizes and complexities.
- Featuring the world's first truly growable contoured cushion, BeMe combines the benefits of contoured seating with the innovative use of patent-pending, high-stretch foam cushion, to maximise pelvic stability and comfort as a child grows and develops. This breakthrough in support technology envelopes the child, increasing pressure distribution and reducing the need for additional accessories such as pommels and hip guides, allowing the child to take centre stage.
- Total ease of use and compatibility are at the heart of BeMe's responsive design. Parts can be easily and intuitively interchanged and clinical adjustments (including seat depth and back height) can be made while the child is seated.
- Armrest, femoral guides and pommels can be quickly released making hoist, standing and side transfers easier than ever.
- Optional air-controlled dynamic backrest which regulates the resistance as a child extends and is speed adjustable as the child returns.



Available in 4 fabric colours

Seat					
	Small	Medium	Large	Medium Long	Large Long
Seat Depth (mm)	170*-325	275-390	350-475	325-440	400-525
Seat Width (mm)	275	325	415	325	415
Max User Weight (kg)**	35	50	70	50	70
Backrest					
	Small	Medium	Large		
Back Height (mm)		360-470	390-520	500-650	
Back Height with Contoured Seat (mm)		370-480	400-530	475-625	
Lateral Width		130-230	200-300	225-350	
Hip Guide Width (mm)		165-230	200-260	250-350	
Backrest Angle		+10 to -25	+10 to -25	+10 to -25	
Max User Weight (kg)**		35	50	70	
Footrest					
Lower Leg length with Basic/Split Seat (mm)		185-285	280-450	290-510	
Lower Leg length with Contoured Seat (mm)		185-285	290-460	300-520	
Maximum Dorsiflexion/Plantarflexion		+/- 10 Degrees	+/- 10 Degrees	+/- 10 Degrees	
Chassis					
Maximum Weight Limit (kg) **		35	50	70	

\*170mm-220mm seat depth required spacer options \*\*System max user weight should be lesser of seat, chassis and backrest



# LECKEY Mygo Seat

The seating system which delivers the optimum in posture, comfort and function for children aged 3-16 years with moderate to complex needs.



## Product Features

- Can be fitted with the newly improved pelvic cradle. The cradle uniquely supports the pelvis by wrapping around the buttocks and posterior pelvis to encourage and maintain an upright functional position. Rigidity has been increased in the back section to enable the cradle to hold its shape, allowing enough flexibility to encourage a posteriorly tilted pelvis into neutral alignment and the overall design prevents users from sliding forward. Extra padding at the belt improves comfort, whilst colour coded buckles and webbing have been added for ease of use.
- Mygo Seat's innovatively designed split seat base enables the legs to move independently beyond the midline to accommodate windsweeping postures- a range of up to 13° abduction and 10° adduction may be achieved. Mygo Seat is also suitable for users with pelvic rotation and tilt, mild scoliosis, shortened hamstrings, femoral and lower leg length difference, and fixed plantar flexion- all this, combined with pelvic stability, ensures alignment of the head and trunk.
- Quickly and easily interfaces with a range of indoor and outdoor hi-low bases, as well as a wide range of Sunrise Medical powered and non-powered mobility bases.



Available in 8 fabric colours and a vinyl version

## Technical Details

	Size 1	Size 2	Size 3
Age [approx]	3-10	8-14	10-16
User Weight	Min 18kg [40lbs] Max 50kg [110lbs]	Min 18kg [40lbs] Max 60kg [132lbs]	Min 18kg [40lbs] Max 70kg [154lbs]
User Height	Min 105cm [41"] Max 150cm [59"]	Min 127cm [50"] Max 168cm [66"]	Min 150cm [59"] Max 175cm [69"]
Seat Width	Min 200mm [8"] Max 325mm [13"]	Min 220mm [8.7"] Max 345mm [13.6"]	Min 250mm [9.8"] Max 400mm [15.7"]
Seat Depth	Min 270mm [10.6"] Max 420mm [16.5"]	Min 350mm [13.8"] Max 470mm [18.5"]	Min 350mm [13.8"] Max 470mm [18.5"]
Chest Width	Min 170mm [6.7"] Max 270mm [10.6"]	Min 170mm [6.7"] Max 270mm [10.6"]	Min 200mm [7.8"] Max 300mm [11.8"]
Backrest Height	Min 360mm [14.2"] Max 470mm [18.5"]	Min 460mm [18.1"] Max 570mm [22.4"]	Min 460mm [18.1"] Max 570mm [22.4"]
Backrest Angle	Prone 10° - Recline 25°	Prone 10° - Recline 25°	Prone 10° - Recline 25°
Lower Leg Length	Min 215mm [8.5"] Max 350mm [13.8"]	Min 315mm [12.4"] Max 470mm [18.5"]	Min 315mm [12.4"] Max 470mm [18.5"]
Armrest Height	Min 160mm [6.3"] Max 210mm [8.3"]	Min 210mm [8.3"] Max 260mm [10.2"]	Min 210mm [8.3"] Max 260mm [10.2"]
Tray Size	550 x 480mm [21.6" x 18.9"]	550 x 480mm [21.6" x 18.9"]	550 x 480mm [21.6" x 18.9"]
Seat Height [footpedal chassis]	340mm-655mm [13.4"-25.8"]	340mm-655mm [13.4"-25.8"]	340mm-655mm [13.4"-25.8"]
Seat Height [powered chassis]	370mm-675mm [14.6"-26.6"]	370mm-675mm [14.6"-26.6"]	370mm-675mm [14.6"-26.6"]



# LECKEY Everyday Activity Seat

Strong, durable, practical and adjustable seating system, providing moderate support for kids and young adults aged 1-18 years.



## Product Features

- Ideal for everyday use- the soft-touch breathable fabric, well-padded ramped seat cushion, contoured back cushion and padded arm rests combine to offer the ultimate in comfort.
- The wide range of accessories and options make it easy to adjust the seat as the user grows and their needs change.
- The flip-up footplate, hi-low base (to floor), long armrests and 10° prone tilt-in-space support self-transfers.
- 4-point pelvic harness gives a secure, stable base of support.
- Height and angle adjustable chest laterals available.
- Easy care – covers can be removed and machine washed.



Available in 3 fabric colours  
and 2 vinyl versions

## Technical Details

	Size 1		Size 2		Size 3	
	mm/kgs	inches/lbs	mm/kgs	inches/lbs	mm/kgs	inches/lbs
Max Load	30kg	66lbs	50kg	110lbs	100kg	220lbs
Weight - Seat & Hi-Low Chassis	29.5kg	65lbs	34kg	75lbs	49kg	108lbs
Seat Height (Hi-Low Chassis)	250-600	9.8-23.6	320-600	12.6-23.6	450-600	17.7-23.6
Seat Width	160-280	6.3-11	260-380	10.2-15	340-500	13.4-19.7
Seat Depth [Standard]	150-280	5.9-11	230-380	9-15	330-480	13-18.9
Seat Depth [Extended]			300-450	11.8-17.7		
Lateral Distance - Chest	150-220	5.9-8.7	170-240	6.7-9.4	240-370	9.4-14.6
Back Height	300-400	11.8-15.7	400-520	15.7-20.5	450-630	17.7-24.8
Back Angle	30° Supine 10° Prone		30° Supine 10° Prone		30° Supine 10° Prone	
Armrest Height	150-250	5.9-9.8	200-300	7.9-11.8	180-325	7-12.7
Lower Leg Length	150-250	5.9-9.8	250-450	9.8-17.7	250-450	9.8-17.7
Lower Leg Angle	-10° to +50°		-10° to +50°		-10° to +50°	
Footplate Angle	-15° to +10°		-15° to +10°		-15° to +10°	
Tilt-in-Space	30° Supine 10° Prone		30° Supine 10° Prone		30° Supine 10° Prone	

# LECKEY Mobility Range

Leckey's interface range enables the BeMe, Squiggles and Mygo Seats to be fitted onto a range of mobility bases from the Sunrise Medical family of manual and power chairs.



## Product Features

- Leckey's interface range enables the Squiggles and Mygo Seats to be fitted onto a range of mobility bases within Sunrise Medical's ZIPPIE and QUICKIE ranges.
- Greater choice and flexibility for Leckey seat users.
- Offers mobility with optimum support and comfort.
- Increased freedom for the user.



Leckey Seat	Sunrise Medical Base	Frame Size Min	Frame Size Max
BeMe Small	Zippie Iris	14"	16"
BeMe Medium	Zippie Iris	17"	18"
BeMe Small	Quickie Iris	14"	16"
BeMe Medium	Quickie Iris	17"	19"
BeMe Large	Quickie Iris	19"	21"
BeMe Small	Quickie SR45	14"	16"
BeMe Medium	Quickie SR45	17"	19"
BeMe Large	Quickie SR45	19"	20"
BeMe Small	Q300	16"	n/a
BeMe Medium	Q400	16"	n/a
BeMe Large	Q400	18"	n/a
Squiggles Seat	ZIPPIE Iris	12"	14"
Squiggles Seat	ZIPPIE TS	12"	14"
Mygo Seat	ZIPPIE Iris	16"	18"
Mygo Seat	ZIPPIE TS	16"	18"
Mygo Seat	QUICKIE SR45	16"	18"
Mygo Seat	QUICKIE TS	16"	18"
Mygo Seat	QUICKIE Salsa M	16"	18"
Mygo Seat	QUICKIE Salsa M2	16"	18"



# LECKEY Leckey Therapy Bench

Developed to support children in a neutral posture, this therapy bench with innovative activity cradles facilitates child active therapy by freeing up the hands.



## Product Features

- Innovative, corseted activity cradle stabilises the pelvis in a neutral position, acting as an extra pair of hands for the therapist whilst giving the child the confidence to move, reach, and grasp.
- Height adjustable trunk support helps the child to sit upright with the head and spine stacked over the pelvis. This is the best and easiest way to learn sitting posture.
- Supportive, height adjustable therapist pad acts as extra support for the therapist who can rest against the rear support and gently encourage shoulder retraction.
- Facilitates lateral and anterior/posterior tilt to promote child activity.
- Height and depth adjustable bench making it versatile for both larger and smaller users.
- Acts as an aid to therapists in the Segmental Assessment of Trunk control (SATCo), providing a greater understanding of where to apply postural support to optimise function.



## Leckey Therapy Bench

Technical Details	
	One Size Only
Max User Weight	80kg (176.3lbs)
Age	1-14
Bench Footprint	620mm x 740mm [63.7" x 29.1"]
Seat Base Height (from ground)	Min 330mm [12.9"] Max 460mm [18.1"]
Seat Base Width	620mm [24.4"]
Seat Base Depth	360mm [14.1"]
Footplate Height (from ground)	Min 20mm [0.7"] Max 110mm [4.3"]
Seat Base to Footplate Height	Min 220mm [8.6"] Max 440mm [17.3"]

## Activity Cradle

Technical Details		
Cradle Size	Hip Circumference (mm)	ASIS - ASIS (mm)
Size 1	Min 400mm [15.7"] Max 500mm [19.6"]	Under 120mm [4.7"]
Size 2	Min 440mm [17.3"] Max 575mm [22.6"]	Min 130mm [5.1"] Max 160mm [6.2"]
Size 3	Min 520mm [20.4"] Max 650mm [25.5"]	Min 160mm [6.2"] Max 180mm [7.0"]
Size 4	Min 600mm [23.6"] Max 740mm [29.1"]	Min 180mm [7"] Max 200mm [7.8"]
Size 5	Min 660mm [25.9"] Max 855mm [33.6"]	Over 200mm [7.8"]

# JCM Star

Multi-functional supportive chair, perfect for use both at home and in the classroom, for children from 1-8 years.

## Product Features

- Can be used directly on the floor for easy access to playtime activities, or with one of the base options for Tilt-in-Space functionality.
- Generous adjustment in the seat length and back height facilitating growth with the child.
- 75mm of height adjustment in each leg allows for growth and also accommodates variable table heights.
- Multi-adjustable arms enable the tray to stay horizontal even when the seat is tilted.
- Built-in pommel helps maintain correct pelvic positioning and can be adjusted in conjunction with the hip pads to align leg position.



Available in black neoprene

Available in 3 Stirata fabric colours

Technical Details		
	Size 1	Size 2
Max User Weight	40kg [88.1lbs]	60kg [132.2lbs]
Seat Depth	Min 200mm [7.8"] Max 300mm [11.8"]	Min 250mm [9.8"] Max 350mm [13.7"]
Seat Height from Floor	Min 285mm [11.2"] Max 360mm [14.1"]	Min 335mm [13.1"] Max 410mm [16.1"]
Seat Width Between Hip Pads	Min 175mm [6.8"] Max 255mm [10.0"]	Min 230mm [9.0"] Max 310mm [12.2"]
Width Between Trunk Lateral Supports	Min 125mm [4.9"] Max 200mm [7.8"]	Min 170mm [6.6"] Max 230mm [9.0"]
Seat to Footplate	Min 130mm [5.1"] Max 270mm [10.6"]	Min 130mm [5.1"] Max 270mm [10.6"]
Footplate Angle	0-30°	0-30°
Armrest Height from Seat	Min 10mm [0.3"] Max 210mm [8.2"]	Min 10mm [0.3"] Max 210mm [8.2"]
Armrest Angle	90° [Prone] 65°	90° [Prone] 65°
Tilt in Space	10° Prone [Star MHT] 15° Tilt [Star MHT]	10° Prone [Star MHT] 15° Tilt [Star MHT]

# JCM Triton

Versatile, modular seating system designed to provide comfort and support for children.

## Product Features

- Numerous options and accessories enable a Triton of any size and specification to be built. The choice of Triton Lite, Triton 0 and Triton 1-3 ensures all ages, from infant to adult, and all postural support needs, from mild to complex, are covered.
- Unique Dynamix options support the user in maintaining an optimum seating position. The unique dynamic back is air pressure, not spring action, controlled, regulating the resistance the user feels as they extend and return in the chair.
- Split seat as standard ensures leg length discrepancy is precisely supported.
- Multi-directional sacral pad to accommodate rotation in the hips.



Available in black neoprene

Available in 4 Libra fabric colours

Available in 3 Stirata fabric colours

Technical Details					
	Triton 0	Triton 1	Triton 2	Triton 2L	Triton 3
Max User Weight	41kg [90.3lbs] [Compact Base] 46kg [101.4lbs] [MHT]	60kg [132.2lbs] [MHT Only]	60kg [132.2lbs] [MHT Only]	80kg [176.3lbs] [MHT Only]	80kg [176.3lbs] [MHT Only]
Seat Depth	Min 150mm [5.9"] Max 275mm [10.8"]	Min 200mm [7.8"] Max 350mm [13.7"]	Min 300mm [11.8"] Max 450mm [17.7"]	Min 415mm [16.3"] Max 565mm [22.2"]	Min 415mm [16.3"] Max 565mm [22.2"]
Seat Width	260mm [10.2"]	315mm [12.4"]	395mm [15.5"]	395mm [15.5"]	475mm [18.7"]
Seat to Footplate	Min 100mm [3.9"] Max 180mm [7.0"]	Min 130mm [5.1"] Max 310mm [12.2"]	Min 200mm [7.8"] Max 380mm [14.9"]	Min 270mm [10.6"] Max 450mm [17.7"]	Min 270mm [10.6"] Max 450mm [17.7"]
Back Angle Recline	5° Prone 40°	5° Prone 40°	5° Prone 40°	5° Prone 25° Recline	5° Prone 40°



# JCM Jupiter

The perfect home seating system for comfort, support and ease of use, catering for children from 3 years right through to adulthood.

### Product Features

- Provides built-in pressure care as standard.
- The combination of Tilt-in-Space and back recline provides maximum comfort by reducing the risk of pressure sores, maintaining posture, increasing circulation to the limbs and easing transfers.
- Height and angle adjustable armrests help to prevent shoulder strain and increase upper body stability.
- The multi-adjustable flip away trunk support, adjustable headrest and leg rests ensure superior comfort and support.



Available in black neoprene

Available in 3 Stirata fabric colours

Available in 4 Libra fabric colours

Technical Details				
	Size 1	Size 2 Standard	Size 2 Deluxe	Size 3
Age	3-9	8-13	8-13	12-Adult
User Weight	60kg [132lbs]	80kg [176lbs]	80kg [176lbs]	100kg [220lbs]
Seat Width	360mm [14.1"]	440mm [17.3"]	440mm [17.3"]	450mm [17.7"]
Seat Depth Single Seat	Min 240mm [9.4"]	Min 360mm [14.1"]	Min 360mm [14.1"]	Min 380mm [14.9"]
	Max 340mm [13.3"]	Max 530mm [20.8"]	Max 530mm [20.8"]	Max 520mm [20.4"]
Tilt-in-Space	-5°	0°	0°	0°
	35°	15°	15°	25°
Back Angle Recline	-5°	-5°	-5°	0°
	25°	35°	35°	25°

Standing  
and Walking



# LECKEY Squiggles+

The Squiggles+ retains the best of the Squiggles Stander but also offers heaps of added benefits so little ones can continue to enjoy standing and be part of the action.



## Product Features

- Anatomically-correct abduction to stabilise the joint and provide optimum loading for early hip joint development. Squiggles+ offers 60° of hip abduction bilaterally in prone, upright or supine. It promotes proper alignment in the hips, knees and feet as pivot points are in line with the hips.
- Increased comfort and support with wrap-around knee supports. The comfortable medial knee supports distribute pressure and are especially suited to abducted standing. The padded wraps maintain the knee in neutral alignment.
- Providing greater confidence with transfers- the larger hip pads improve the carer's and child's confidence and feeling of security when transferring. The improved pivot chassis facilitates full horizontal loading.
- Activity tray with improved height and adjustability.
- Height-adjustable, foldable push handle.



Available in 3 fabric colours and a vinyl version



Technical Details	
One Size Only	
Age	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5 inches Max 1110mm / 43.5 inches
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches
Chest Support Height (footplate to armpit)	Min 460mm/18.1 inches Max 820mm / 32 inches
Chest Width (distance between laterals)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Hip Support (height from footplate to mid buttocks)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Hip Width (distance between hip guides)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Distance between Midline of Knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°



# LECKEY Squiggles TT

Evolutionary early intervention standing frame that facilitates Targeted Training for children aged 1-5 years.



## Product Features

- Maintaining the great features of the Squiggles+, the Squiggles TT goes a step further and facilitates Targeted Training. This evidence-based therapy is carried out with the child in an upright position, training one segment of their trunk at a time, allowing them to gain sequential control from the top down. Improved trunk control will help increase functional skills e.g. self-care and play. Better head control can help stop a child dribbling, improve their eye control as well as their ability to communicate.
- Squiggles TT is the only product specifically designed for Targeted Training therapy and which caters for all children. Horizontal loading is possible for complex children with poor head control and Squiggles TT is the only product to target all 7 trunk segments.
- Offers a fun and easy way for parents to carry out therapy at home. The snug supportive wraps are easy to use, the ratchet system with repeatable locking eliminates the fear of changing settings and the play ideas booklet provides support for home use.
- Dynamic base allows children to experience rocking movement whilst gaining upright postural control.
- Padded tray facilitates Targeted Training at the upper thoracic level.
- Conversion kit available for existing Squiggles+ users.



## Technical Details

	One Size Only
Age (approx.)	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5 inches Max 1100mm / 43.3 inches
Trunk Support Height (Distance from footplate to armpit)	Min 160mm / 6.25 inches Max 820mm / 32 inches
Trunk / Pelvic Width	Min 135mm / 6.25 inches Max 250mm / 9 inches
Pelvic Support Height (Distance from footplate to mid-buttock)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Distance between midline of knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Hip Width (distance between hip guides)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Distance between Midline of Knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate	Plantarflexion -10° Dorsiflexion -10°
Pivot chassis angle range	90-180°
Tray Height	Min 480mm / 18.5 inches Max 825mm / 32.5 inches



Available in 3 fabric colours and a vinyl version

# LECKEY Squiggles Stander

Versatile, supportive, compact stander which offers the perfect start to standing for little ones.



## Product Features

- 3-in-1 prone, supine and upright stander.
- Ideal for home use due to ease of adjustments.
- Minimalist child-focused design with fun bright pads blends easily into nursery/school environments.
- Low to the ground, encouraging socialisation and interaction at eye level.
- Individually adjustable chest and pelvic supports are cushioned for comfort and adjust easily in height, width and angle to suit individual needs, ensuring a comfortable and secure position.
- Frame can be easily disassembled for storage and transportation.



Available in 3 fabric colours

Technical Details	
	One Size Only
Age	1-5
User Height	730mm-1100mm [28.5"-43.3"]
Max. User Weight	22kg [48lbs]
Chest Support Height (footplate to armpit)	550mm-820mm [21.6"-32.2"]
Chest Width (distance between laterals)	160mm-230mm [6.2"-9"]
Chest Support Adjustment	+/- 30° 50mm [1.9"] - Depth
Hip Support (height from footplate to midbuttocks)	300mm-530mm [11.8"-20.8"]
Hip Width (distance between hip guides)	160mm-230mm [6.2"-9"]
Distance between Midline of Knees	140mm-215mm [5.5"-8.4"]
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°
Tray Height (footplate to top of tray)	555mm-825mm [21.8"-32.4"] - Prone 480mm-825mm [18.8"-32.4"] - Supine
Tray Angle Adjustment	30° - Prone 40° - Supine



# LECKEY Mygo Stander

The stander which offers infinite positioning, helping children with moderate to complex needs stand with confidence.



### Product Features

- 3-in-1 prone, supine and upright stander.
- Designed to improve upper limb function. The chest pad design replicates the shape of the rib cage which allows the shoulder girdle to freely protract and retract in prone, whilst in supine, the shoulder support promotes shoulder protraction. Furthermore, Mygo offers a lower tray height than most standers.
- Accommodates hip and knee flexion contractures of up to 25 degrees. The movement of the footplate and knee supports are linked which allows the user to stand with hips and knees flexed. This results in a more upright posture combined with the ability to support users whose contractures previously made them unsuitable for standing therapy.
- Use of the pelvic derotation belt puts the user's pelvis in a central and neutral position with minimum effort.



Available in 4 fabric colours

Technical Details		
	Size 1	Size 2
Age	4-10	8-14
User Height	1080mm-1450mm (42.5"-57")	1350mm-1700mm (53.1"-66.9")
Max. User Weight	50kg (110lbs)	60kg (132lbs)
Chest Support Height (footplate to armpit)	760mm-1080mm (29.9"-42.5")	980mm-1300mm (38.5"-51.1")
Chest Width (distance between laterals)	160mm-280mm (6.2"-11")	200mm-320mm (7.8"-12.5")
Hip Support Height (footplate to mid-buttocks)	460mm-730mm (18.1"-28.7")	620mm-875mm (24.4"-34.4")
Hip Width (distance between hip guides)	220mm-330mm (8.6"-12.9")	280mm-390mm (11"-15.3")
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°	Plantarflexion 10° Dorsiflexion 10°
Max. Knee Contracture Accommodated	25°	25°
Tray Height Range (from top of hip pad)	0mm-150mm (0"-5.9")	0mm-150mm (0"-5.9")
Tray Angle Adjustment	45° - Prone 45° - Supine	45° - Prone 45° - Supine



# LECKEY Horizon

Robust, durable 3-in-1 stander for users aged 4 years to adult.



### Product Features

- Accommodation of leg length discrepancy, scoliosis, lumbar lordosis, kyphosis, knee flexion contractures and asymmetrical head position. The frame's growth potential and the wide range of accessories accommodate changing need and potential asymmetry as the user grows.
- Perfectly sized for hoisting and supine standing. The large landing space for hoisted transfers coupled with support pads which provide a proximal fit to suit each user provides the perfect combination.
- Robust, stable and durable frame ensures ongoing confident product use.
- Individually adjustable chest, hip and knee pads as well as footplates provide unsurpassed levels of adjustability.
- Height adjustable hip laterals and pelvic belt, enabling de-rotation of the hips, alignment of the pelvis and centre of gravity over the base of support.



Available in 3 fabric colours

Technical Details			
	Size 1	Size 2	Size 3
Age	4-10	8-14	14-adult
User Height	1000mm-1400mm [39.3"-55.1"]	1400mm-1800mm [55.1"-70.8"]	1550mm-1800mm [61"-70.8"]
Max. User Weight	50kg (110lbs)	80kg (176lbs)	100kg (220lbs)
Shoulder Height (footplate to top of pad)	820mm-1150mm [32.2"-45.2"]	1000mm-1450mm [39.3"-57"]	1000mm-1450mm [39.3"-57"]
Hip Pad Height (footplate to top of pad)	450mm-750mm [17.7"-29.5"]	750mm- 950mm [29.5"-37.4"]	750mm-1150mm [29.5"-45.2"]
Hip Pad Width	175mm-270mm [6.8"-10.6"]	210mm-390mm [8.2"-15.3"]	310mm-470mm [12.2"-18.5"]
Knee Height (footplate to middle of pad)	250mm- 450mm [9.8"-17.7"]	400mm-550mm [15.7"-21.6"]	400mm-650mm [15.7"-25.5"]
Width Between Chest Laterals	200mm-300mm [7.8"-11.8"]	240mm-380mm [9.4"-14.9"]	340mm-460mm [13.3"-18.1"]
Footplate Angle Adjustment +/-	10°	10°	10°
Angle Adjustment	0°-90°- Prone 0°-90°- Supine	0°-90°- Prone 0°-90°- Supine	0°-90°- Prone 0°-90°- Supine
Tray Angle Adjustment	0°-40°- Prone +/- 20°/30°- Supine	0°-40°- Prone +/- 20°/30°- Supine	0°-40°- Prone +/- 20°/30°- Supine



# LECKEY Pronestander

Robust, adjustable prone standing system for users aged 1-18 years.

## Product Features

- Flip-up sandals support users completing a standing transfer directly into the product.
- Long base ensures maximum stability for even the most active users, providing caregivers with peace of mind.
- Unique de-rotation belt and pelvic pad work together to keep the user in a central position within the frame with equal pressure over each foot.
- Individually adjustable chest, hip, knee and sandal supports to meet a wide range of clinical needs.
- Easy to set up and adjust with five incremental secure prone positions from 85 to 45 degrees.



Size 1 only available in Meadow Fabric  
Size 2 & 3 only available in Mondeo Fabric

Technical Details			
	Size 1	Size 2	Size 3
Age	1-5	4-10	9-18
User Standing Height	800mm-1151mm [31.4"-45.3"]	1100mm-1448mm [43.3"-57"]	1400mm-1849mm [55.1"-72.7"]
Max. User Weight	60kg [132lbs]	80kg [176lbs]	90kg [198lbs]
Chest Pad Height	470mm-680mm [18.5"-26.7"]	680mm-975mm [26.7"-38.3"]	975mm-1350mm [38.3"-53.1"]
Prone Angles	85°, 75°, 65°, 55°, 45°		
Width between Chest Laterals	152mm-254mm [5.9"-10"]	203mm-305mm [7.9"-12"]	254mm-356mm [10"-14"]
Tray Sizes	540mm x 385mm [21.2" x 15.1"]	580mm x 410mm [22.8" x 16.1"]	645mm x 460mm [25.3" x 18.1"]

# LECKEY Freestander

Robust, upright standing system for users aged 1-18 years.

## Product Features

- Flexibility to alter the frame to meet individual user needs e.g. pelvic support can be removed for users with trunk control, conversely, the tray height and angle can be increased where more support is required. A simple, active set-up is also possible for more able users e.g. hemiplegic or GMFCS II/III.
- Straightforward, intuitive adjustment as all supports run off two vertical uprights.
- Pelvic belt achieves a unique position which enables de-rotation of the hips, alignment of the pelvis and centre of gravity over the base of support as well as accommodation or correction of flexion at the hips and knees.



Size 1 only available in Meadow Fabric. Size 2 & 3 only available in Mondeo Fabric.

Technical Details			
	Size 1	Size 2	Size 3
Age	1-5	4-10	9-18
User Standing Height	800mm-1151mm [31.4"-45.3"]	1100mm-1448mm [43.3"-57"]	1400mm-1849mm [55.1"-72.7"]
Max. User Weight	60kg [132lbs]	80kg [176lbs]	90kg [198lbs]
Chest Pad Height	470mm-680mm [18.5"-26.7"]	680mm-975mm [26.7"-38.3"]	975mm-1350mm [38.3"-53.1"]
Width between Chest Laterals	152mm-254mm [5.9"-10"]	203mm-305mm [7.9"-12"]	254mm-356mm [10"-14"]
Tray Sizes	540mm x 385mm [21.2" x 15.1"]	580mm x 410mm [22.8" x 16.1"]	645mm x 460mm [25.3" x 18.1"]

# LECKEY Totstander

Simple and easy to use, upright standing system for kids aged 1-5 years.

### Product Features

- Can be used as a stander or by therapists for standing assessments.
- Promotes active standing, movement and weight shift.
- Compact and easily transported to different environments e.g. from therapy centres to the home.
- Provides a gentle introduction to rehab equipment via a simple, easy-to-use, child-friendly design.



Available in 1 fabric colour combination

Technical Details	
	One Size Only
Age	1-5
Chest Band Height	470mm- 640mm (18.5"-25.1")
Base Depth	450mm (17.7")
Base Width	450mm (17.7")
Max. User Weight	22kg (48lbs)

# Timoteos Ella

4-in-1 stander which accommodates a full range of postural positions, including abduction.

### Product Features

- Accommodates users from 600mm to 1800mm across 4 sizes. An optimal fit can be achieved for users through interchangeable sizes of accessories, stepless adjustments, asymmetrical adjustments and ability to fine tune set-up whilst the client is in the stander.
- Only one tool is needed for setup and adjustments- adjustments are intuitive and easy to do, supporting a user-friendly approach to postural management.
- The comfort offered by Ella is second to none soft padding and comfortable upholstery mean you can be confident the user is always experiencing great comfort.
- Individual abduction is available for each leg. Furthermore, the abduction brackets can be adjusted independently, meaning asymmetric posture can be accommodated.
- Available with a choice of manual or power base which can be detached for easy transportation and storage.



Available in 3 Stirata fabric colours

Available in 4 Dartex wipeable fabric colours

Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	0-4	3-7	6-12	11-18
User Height	600-950mm	900-1300mm	1200mm-1550mm	1300mm-1800mm
Chassis Footprint Width	580mm	580mm	702mm	752mm
Chassis Footprint Length	750mm	750mm	1000mm	1180mm
Trunk Support Pad (Top) to Footplate Height	430-680mm	570-875mm	800-1280mm	900-1500mm
Pelvic Support Pad (Top) to Footplate Height	300-540mm	475-785mm	500-1040mm	625-1225mm
Knee Cup (Top) to Footplate Height Angle Adjustment	140-450mm	180-550mm	200-650mm	250-850mm
Maximum User Weight	30kg	40kg	55kg	80kg



# LECKEY MyWay+

MyWay+ is an open frame walker which cleverly supports children in an upright position, improving stepping and maximising opportunities for exploration, interaction and participation.



## Product Features

### • Tailored Approach To Every Child

With the wide range of adjustments and accessories, the MyWay+ can be tailored to benefit nearly all children [age 1-16] who need additional trunk support. The robust new multi-functional castors, with indoor/outdoor performance will open up a world of possibilities.

### • Unique Corset Like Harness

The perfectly engineered, corset like harness provides 360° wraparound support to the trunk. This sensory feedback supports low tone, dampens high tone and ultimately improves stepping pattern. Hoisting loops make transferring in and out a breeze.

### • Open Frame Design

Promoting freedom and independence are at the heart of MyWay+. It's what it's designed for. The open-frame design ensures easy, peer-to-peer interaction and increases participation.



S1 S2 S3 S4 S5 S6

Harness available in 6 fabric colours

	Size 1 Frame	Size 2 Frame	Size 3 Frame
Age [approx.]	1-5	4-11	9-16
User Shoulder Height [to floor] [mm]	590 – 870	870-1190	1070 – 1490
User Saddle Height [inside leg] [mm]	240 – 410	400-610	550 - 810
Max. User Weight [kg]	25	50	80
Frame Weight [kg]	14 [S1 frame and upper support unit do not separate]	7.5	11
Upper Support Unit Weight [kg]		10	11
Frame Width [mm]	570	670	720
Frame Length [mm]	840	930	1180

Harness Size	Age (Approx)	Waist at Navel Level [cm]	Max. User Weight [kg]	Colour
1	1-4	46-56	50kg	Purple
2	2-6	50-60	50kg	Green
3	5-9	55-65	50kg	Blue
4	8-12	60-70	80kg	Red
5	12-15	70-80	80kg	Orange
6	16-18	80-90	80kg	Grey

# LECKEY MyWay Pedal

MyWay Pedal is a revolutionary elliptical trainer which helps improve strength and functional ability. Attaching to the size 2 and 3 MyWay frame, MyWay Pedal is ideal for giving on-the-go kids the opportunity to be active, have fun and join in with their friends.



## Product Features

- Targets and strengthens anti-gravity muscles to help reduce the risk of hip and knee contractures developing and improve function. Pedalling in an upright extended position targets and strengthens anti-gravity muscles (calf, quads and hamstring)-unlike a trike which encourages a flexed position, the upright position of MyWay Pedal ensures the quads are worked exactly where strength is often lacking, at the end of range. Working the end of range will reduce the risk of hip and knee contractures and limit the potential for deterioration in standing and walking ability.
- Facilitates elliptical action more like stepping to help support gait therapy goals. MyWay Pedal provides a sensori-motor experience similar to gait. Unlike the circular motion of cycling, the elliptical pedalling of MyWay Pedal gives an extended stretch of the lower limbs as it involves a longer forwards/backwards movement rather than an up/down movement. It helps improve loading around the hip to aid joint development.



- Use in conjunction with the Leckey Connect App to chart performance and motivate the child. The app provides a valuable way of measuring performance and progress on a child's speed, distance and duration of elliptical pedalling. It is available on Android Smartphone, tablet and the Leckey online portal.
- Attaches to the size 2 and 3 MyWay frame and can be fitted with the user in or out of the frame.
- Two sandal sizes available.
- Height adjustable and removable handlebar.





# Say hello to Firefly



# firefly GottaGo

Toileting Made Easy, Home and Away

An innovative, portable toilet seat which makes toilet training easier for families of children with special needs.



## Product Features

- Suitable for children aged 2-9 years with a maximum user weight of 30kg/66lbs.
- GottaGo frame has 3 height settings: potty height, standard toilet height and disabled toilet height.
- Warm-touch, wipe clean materials guarantee a comfortable and snug fit.
- Easy to set up and put away.
- Specially tailored toilet training guide to support and motivate families throughout their toilet and potty training journey.



## Technical Details

	Size 1	Size 2
Age (approx.)	2-5 years	4-9 years
Hip Width	22.5cm	27cm
Minimum Seat Depth*	22cm	29cm
Backrest incline angle (forward)	8°	25°
Thigh angle (off horizontal)	40°	40°
Backrest height	24cm	36cm
Weight (seat and frame)	5.8kg	6.1kg

\* For boys sizing it is recommended to choose the larger size providing they meet the minimum seat depth measurement.



# firefly GoTo Seat

One Little Seat, Countless Possibilities

A lightweight, portable and stylish postural support seat that enhances children's lives through participation.



## Product Features

- Suitable for attachment to most standard chairs in the home, restaurants and classroom by using the adjustable attachment straps. Can also be attached to shopping trolleys.
- Adjustable headrest and laterals that adjust for growth.
- Standard or Advanced Headrest available.
- Pressure relief cushion available.
- Optional floorsitter accessory has 26 recline positions, from upright (106°) to right back (140°).



Available in 2 fabric colours and a vinyl version

## Technical Details

	Size 1	Size 2
Age	1-3	3-8
Seat Base to Top of Headrest	380mm-540mm (14"-21")	510mm-680mm (20"-27")
Lateral Support Width	190mm-260mm (7.5"-10")	210mm-300mm (8"-12")
Seat Depth	175mm (7")	244mm (9.5")
Max. User Weight	15kg (33lbs)*	30kg (66lbs)**

\*If a child is above 13kg/28.5lbs GoTo Seat cannot be used in a shopping trolley.

\*\*Size 2 GoTo Seats are too large for use in shopping trolleys.



# firefly Upsee

Little Steps to an Unexplored World

A mobility harness designed to allow children with special needs to experience upright movement like anyone else.



## Product Features

- Facilitates gradual weight-bearing and gives children with motor impairment the closeness and confidence needed to stand and move with the help of an adult.
- Creates opportunity for social and emotional development as users can interact with peers at eye-level.
- Comprised of three parts - child harness, adult belt and double sandal.
- Durable, machine washable and breathable fabrics.
- Two colour options available for each size - please visit our website for further details.



Available in 4 fabric colours

## Technical Details

	Extra Small	Small	Medium	Large
Age	1-2	2-4	4-6	6-8
Max. User Weight	15kg (33lbs)	15kg (33lbs)	20kg (44lbs)	25kg (55lbs)
Shoulders to Hip	280mm-320mm (11"-13.5")	310mm-360mm (12"-14")	350mm-390mm (14"-15.5")	380mm-420mm (15"-16.5")
Chest Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	570mm-650mm (22"-26")	630mm-710mm (25"-28")
Lower Hip Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	580mm-670mm (23"-26")	640mm-730mm (25"-29")



# firefly Scooot

## Stylish Wheels for Little Explorers

A 4-in-1 sit and ride which gives children the freedom to explore and have fun in style.



### Product Features

- Scooot has four configurations: Scoot, Ride, Pooosh and Crawl. Different configurations for different abilities and different ways to play.
- For added support in Scooot, the Advanced Backrest Cushion helps. It's an optional extra that gives support to the child's trunk. The laterals move up and down to give support where it's needed. Or remove them and just have a comfy cushioned backrest. The Advanced Backrest Cushion is easy to attach, wash, and can be bought separately.
- One size caters for children aged 2-6- the footrest adjusts to suit the child's leg-length as they grow.



### Four Configurations

- **Crawl** allows kids to move and explore on their tummy, helping to strengthen their upper body and tummy muscles as they lift their heads and shoulders.
- **Scoot** lets kids sit upright to move around using their feet. Supported by a lap belt and backrest, the child's trunk and pelvis are encouraged to stabilise, freeing up their hands to lift and handle objects, play with family, and explore their surroundings.
- **Ride** lets kids propel themselves around the floor by moving the wheels with their hands. Having fun in this function can build upper body strength, whilst exploration and participation can improve spatial awareness as well as cognitive, sensory and communication abilities.
- **Pooosh** lets somebody else do the hard work. If the child isn't ready to propel themselves yet (or they're just tired), attach the Pooosh handle. It helps get little ones used to the Scooot. And it means parents and therapists can help them follow the fun at home, school or therapy sessions.

### Technical Details

Size	4-in-1
Age	2-6
Max. User Weight	22kg (48.5lbs)
Seat Depth	400mm (15.5")
Seat Width	230mm (9")
Hip to end of footplate	520mm-630mm (20.5"-24.5")
Backrest Width	260mm (10.5")
Backrest Height	300mm (12")
Wheel Arch Width	230mm (9")
Pooosh Pole	800mm-890mm



# firefly Playpak

Fun Therapy for Little Learners

A portable, versatile early intervention therapy kit which provides children with the building blocks they need to help develop key skills.



## Product Features

- Lightweight and portable.
- More than just practical- Playpak looks great with its colourful, fun, child-friendly design.
- Removable fabric can be washed at 40°C and tumble dried at a low temperature.
- Online Therapy Booklet, featuring expert advice and instructions to help develop physical, cognitive, sensory and communication skills.
- Lightweight bag containing an entire activity centre, including rolls, wedges and supports (designed by clinical experts) that stick together in dozens of different ways to practice the early positions which form the building blocks for later abilities.



## Helps develop 5 key skills:

- Back Lying (Supine)
- Tummy Lying (Prone)
- Floor Sitting (Long Sitting)
- Side Lying
- Four Point Kneeling





# firefly WeGo

Everyday Adventures in Style

A stylish, lightweight pushchair which is perfect for everyday use by active families and little ones who like to be on the move.



## Product Features

- Compatible with the Firefly GoTo Seat.\*
- Lightweight and easy to fold.
- Accessories included as standard- sun canopy, foot muff, rain cover and basket.
- Recline angle up to 140 degrees.
- Maximum user weight of 30kgs.\*
- Wheels suitable for rough terrain.
- Small footprint so can fit in most car boots.

\* Please check sizing guide for exact measurements



## Technical Details

Footprint Length (Collapsed)	116cm
Footprint Height (Collapsed)	38cm
WeGo Height	122cm
Seat Depth	26cm
Seat Height	50cm
Seat Width	35cm
Footrest Shortest Setting	25cm
Footrest Longest Setting	33cm
WeGo Weight	11kg
WeGo Weight with GoTo Seat	13kg

Please note - Not all children using the GoTo Seat Size 2 will be suitable for the WeGo. Please refer to our sizing guide.

# firefly Splashy

## Water Adventures for Little Splashers

A lightweight, portable multipurpose bath seat that can be used for both bathing and other fun water play activities.



### Product Features

- Splashy is designed to sit low in the bath to allow the child to benefit from being immersed in the water.
- Great aesthetics via bright colour options and a modern design.
- Splashy adapts to the child's postural support needs, with the option of a 5-point harness or the 3-point pelvic support. There are four detachable bumpers to add extra support wherever it's needed. The harness and foam supports attach easily to the seat base using toggles.
- Warm-touch foam feels soft whether it's wet or dry.
- Splashy's base unit lets you switch easily between 26 recline positions, from upright (106 degrees) to full recline (140 degrees).
- Lightweight, detachable parts make it easy to clean.



Available in 3 colours

### Technical Details

Age	1-8
Backrest Height	610mm (24")
Backrest Width	385mm (15")
Backrest Depth	300mm (12")
Seat Width	350mm (14")
Floorsitter Width	420mm (16.5")
Floorsitter Length	630mm (24")
Max Weight	30kg (66lbs)



firefly

Splashy Big

Water Adventures for Big Splashers

With all of the features and benefits of Splashy, the Splashy Big opens up many more opportunities for 8-16 year olds who want to participate in fun activities with friends and family.

Product Features

- Suitable for 8-16 years with a max user weight of 55KG/122lbs.
- Splashy Big is lightweight and super easy to maintain with its warm touch and quick drying covers. Cleaning it down after use is a breeze.
- Splashy Big's floorsitter has 26 recline options plus a new stabiliser feature for use outside of the bath.



Available in 2 colour combinations



Technical Details	
Seat Backrest Height	720mm [28.5"]
Seat Backrest Width	450mm [17.5"]
Seat Base Depth	400mm [15.5"]
Seat Base Width	450mm [17.5"]
Floorsitter Width	420mm [16.5"]
Floorsitter Width with Stabiliser	665mm [26"]
Floorsitter Length	795mm [31"]
Maximum Weight	55kg [121lbs]
Age Range	8-16 years old



Other  
Products

# LECKEY Advance Bath Chair

The Leckey Advance Bath Chair provides postural support while ensuring the individual is safe and secure during bathing.

## Product Features

- Height adjustable legs.
- Breathable cover, removable and machine washable.
- Back angle adjustment from 0-90 degrees.
- Folds flat for storage



Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	1-5	4-9	8-14	12-18
Min User Height	750mm (30")	950mm (37")	1150mm (45")	1350mm (53")
Max User Height	1050mm (41")	1250mm (49")	1450mm (57")	1650mm (65")
Max User Weight	72kg (159lbs)	72kg (159lbs)	72kg (159lbs)	72kg (159lbs)



Available in 2 fabric colours

# LECKEY Shower Trolley

The Shower Trolley complements the Advance Bath Chair which can be converted within seconds into a shower chair.

## Product Features

- Converts the Advance Bath Chair to a versatile shower chair in seconds.
- Manufactured from medical grade stainless steel.
- Robust frame and high quality lockable castors.
- Can be easily and safely manoeuvred over wet and slippery surfaces.
- Can be adjusted to fit all four sizes of the Leckey Advance Bath Chair.



Technical Details	
	One Size Only
Age	1-18
Length	900mm-1150mm (36"-45")
Width	600mm (24")
Seat Height from Floor	650mm-700mm (26"-28")

# LECKEY Sleepform

Strong, durable, practical and adjustable sleeping system for kids aged 0-18 years.

## Product Features

- Close proximal fit cocoons the child and makes them feel safe and secure. Adapts quickly and easily to changing body shape.
- Compact and portable for visits to family, going on holiday etc. Suitable for use on a variety of surfaces.
- Supports 24-hour postural management as can be used during the day to support changes of position.
- Can be used in cots and single beds and suits all sleeping positions.
- Air flow mattress.
- Grid sheet to record best sleeping positioning.
- Temperature control sheet available.



Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	0-1	1-5	1-5	5-18
Suitable for	Cot	Cot	Single Bed	Single Bed



# WHITMYER® Headrests

The Whitmyer range provides a variety of head positioning solutions for support and/or control of the cervical spine and head, from the most basic of needs, to the most complex. The range offers comprehensive sizing and accessories, allowing for individualised positioning of the user.

## Adjust-A-Plush and Plush

### Product Features

- Single pad headrests which are easily adjustable and accommodate changes or inconsistencies in body positioning.
- Utilises dual density foam layers on a contoured, curved panel.
- Unique “cupping” shape provides posterior support with a high degree of comfort and safety.
- Basic headrest which is ideal for tilt/ recline and has the ability to add lateral and anterior support. The larger size pads are ideal for users with active movements who may try to hook their heads behind other headrests.



Technical Details					
	Adjust-A-Plush Standard	Adjust-A-Plush Narrow	Plush 250mm	Plush 350mm	Plush 480mm
Width	200mm-460mm [7.8"-18.1"]	152mm-355mm [5.9"-13.9"]	235mm [9.2"]	336mm [13.2"]	182mm [18.9"]
Height	127mm [5"]	76mm [2.9"]	127mm [5"]	127mm [5"]	190mm [7.4"]

## Contoured Cradle

### Product Features

- Unique anatomic shape, designed to cradle the occiput and provide lateral cervical support, supplying a high degree of comfort and positioning.
- Single rotational collar adjustment allows positioning for maximum support surface contact.
- Multiple pad sizes fit infant to adult



Technical Details			
	Small	Medium	Large
Width	216mm [8.5"]	250mm [9.8"]	330mm [12.9"]
Height	100mm [3.9"]	120mm [4.7"]	176mm [6.9"]

## Heads Up

### Product Features

- Innovative anterior head support system.
- Designed to maximise users' interaction with their environment by providing anterior, posterior and lateral support of the head and neck.
- Features a chin prompt as well as swing-out, adjustable arms to easily accommodate transfers and varying user presentations.



Technical Details								
	Early Intervention		Small		Medium		Large	
	Short	Long/Tall	Short	Long/Tall	Short	Long/Tall	Short	Long/Tall
Width	84mm [3.3"]	84mm [3.3"]	102mm [4"]	102mm [4"]	117mm [4.6"]	117mm [4.6"]	147mm [5.7"]	147mm [5.7"]
Height	64mm [2.5"]	91mm [3.5"]	71mm [2.7"]	109mm [4.2"]	84mm [3.3"]	127mm [5"]	94mm [3.7"]	150mm [5.9"]

## Notes



LECKY



ALLIED  
MEDICAL

[helpis@alliedmedical.co.nz](mailto:helpis@alliedmedical.co.nz) | [www.alliedmedical.co.nz](http://www.alliedmedical.co.nz) | 0800 31 61 81